

**LUDLUM MODEL 43-10-6
ALPHA-BETA SAMPLE COUNTER**

**SERIAL NUMBER PR151543
AND SUCCEEDING SERIAL NUMBERS**

December 2015

**LUDLUM MODEL 43-10-6
ALPHA-BETA SAMPLE COUNTER**

**SERIAL NUMBER PR151543
AND SUCCEEDING SERIAL NUMBERS**

December 2015



LUDLUM MEASUREMENTS, INC
501 OAK STREET, P.O. BOX 810
SWEETWATER, TEXAS 79556
325-235-5494, FAX: 325-235-4672

**Model 43-10-6 Alpha-Beta Sample Counter
December 2015**

TABLE OF CONTENTS

1. GENERAL	1
2. OPERATION	1
3. SPECIFICATIONS	2
4. CALIBRATION PROCEDURE	2
4.1 Counting Instrument	2
4.2 Procedure	2
5. TROUBLE SHOOTING	4
5.1 Zero or Very Low Counts	4
5.2 No Source Plateau	4
5.3 Excessive Background Count	4
6. PARTS LIST	5
Model 43-10-6 Alpha Sample Counter	5
Assembly View, Drawing 142 x 113	5
Switch Filter Board, Drawing 142 x 58	5
Voltage Divider Board, Drawing 435 x 964	5
DRAWINGS AND DIAGRAMS	6

Model 43-10-6 Alpha-Beta Sample Counter

December 2015

1. GENERAL

The Model 43-10-6 is an Alpha-Beta Sample Counter capable of holding up to a 7 cm (2.75 in.) diameter filter. Use of the filler plate allows smaller samples to fit snugly in the sample drawer.

The sample drawer, when fully closed, strikes a microswitch to allow high voltage to be applied to the photomultiplier tube (PMT). The sample drawer is locked in the closed position by rotation of the slide lever mounted on the side of the instrument. To discriminate between alpha and beta radiation simultaneously, the counting instrument has

to have either separate power supplies or threshold controls for each channel. The Ludlum Model 2929 Scaler, Model 2223, or Model 2224 instruments provide the necessary circuitry for simultaneous alpha-beta discrimination. ZnS (Ag) is used for alpha radiation detection, and plastic scintillation material is used for detection of beta radiation. The scintillation material is covered by 0.4 mg/cm² metallized polyester to reduce light response (excessive background).

2. OPERATION

- Connect the Model 43-10-6 to the scaler. The coax cable with "C" connectors carries both the signal and high voltage (HV).
- HV is applied to the photomultiplier tube (PMT) when the sample drawer is pushed completely in, tripping the microswitch. Rotate sample drawer lever to the locked position, securing sample drawer in the ON position.
- Alpha background count is approximately 3 counts per minute (cpm) or less. Beta background count is approximately 80 cpm.
- To count a radioactive sample, place sample on appropriate side of sample holder for the one- or two-inch filters. Do not allow sample to extend above the top of the sample slide.
- A background count should be taken after each source count to check for contamination on sample holder or area within O-ring.

Model 43-10-6 Alpha-Beta Sample Counter
December 2015

3. SPECIFICATIONS

Scaler Input Sensitivity: 2 to 200 mV

Photomultiplier Tube: 5.1 cm (2 in.) diameter, 14-pin tube base, 10-pin dynode structure

Scintillator Material: ZnS disc; plastic (EJ212) 0.025 cm (0.01 in.) thick

Detector Operating Voltage: 500-1200 Vdc

Window: 0.4 mg/cm²

HV Switch: opening sample slide disables PMT high voltage

Sample Sizes: 1 1/8 inch diameter x 3/16 inch deep, 2 inch diameter x 3/16 inch deep, 2-1/4 inch diameter x 3/8 inch deep, or 2-3/4 inch diameter X 1/8 inch deep

Efficiency (4 pi): 35% for ²³⁹Pu, 5% for ¹⁴C

Channel Cross Talk: alpha in beta channel ≤10%; beta in alpha channel ≤1%

Construction: aluminum housing with beige powder cotang.

Size: 23.6 x 11.4 x 23.6 cm (9.3 x 4.5 x 9.3 in.) (H x W x L)

4. CALIBRATION PROCEDURE

Caution: Do not tip sample counter over with sample holder in sample slide. The sample holder will tear the thin metallized polyester window, allowing light to scintillate the ZnS, which causes excessive counts in the beta channel.

For instruments with separate power supplies (fixed threshold), the alpha channel will operate at a lower voltage than the beta channel.

4.1 Counting Instrument

For Calibrated Scaler Instrument:

HV Range is Nominally 800 ±200 volts. For instruments with separate power supplies and fixed threshold, the alpha channel will operate at a lower voltage than the beta channel.

Nominal input sensitivity and alpha channel at 175 mV, beta channel at 4 mV, and upper discriminator set at 50 mV

The Model 43-10-6 instrument should be dark-adapted before taking data.

4.2 Procedure

Determining Operating Voltage

- Connect the Model 43-10-6 to the counting instrument with proper cable.
- Place a calibrated ¹⁴C source in the sample holder. Close and lock the sample drawer.
- Adjust counting instrument HV until it receives 5% (or greater) 4pi efficiency is achieved.
- Decrease HV by 25 volts.
- Record the HV.
- Record ¹⁴C source count and beta crosstalk in alpha channel.

Model 43-10-6 Alpha-Beta Sample Counter
December 2015

- Remove ^{14}C source and record background count in alpha and beta channels.
- Place calibrated ^{239}Pu source in sample holder. Close and lock sample drawer.
- Record ^{239}Pu source count and alpha crosstalk in beta channel.
- Increase the HV by 25 volts.
- Repeat the six previous steps until one or more of the following conditions are met:
 1. Beta background exceeds 80 cpm.
 2. Alpha background exceeds 3 cpm.
 3. Alpha crosstalk in beta channel exceeds 10%.
 4. Beta crosstalk in alpha channel exceeds 1%.
- The operating voltage should be selected as a point where:
 1. ^{14}C efficiency $\geq 5\%$
 2. ^{239}Pu efficiency $\geq 37\%$
 3. crosstalk alpha in beta channel $\leq 10\%$
 4. crosstalk beta in alpha channel $\leq 1\%$

Determining Efficiency

- NIST traceable check sources required.
- Set HV as determined previously.
- Record a one-minute background and one-minute source count. Subtract the background from the source count. Divide the net source count by the dpm value of the source, times 100 for 4π efficiency.

If the source value is listed in microcuries (activity):

- Convert the microcurie value to a dpm value by multiplying the microcurie value by 2.22×10^6 . Calculate the 4π efficiency as in the previous steps.

5. TROUBLE SHOOTING

5.1 Zero or Very Low Counts

- Large light leak
- PMT malfunction
- Broken wire in tube socket
- Inoperative HV switch on Model 43-10-6 or broken wire
- Counting instrument malfunction
- Source too far from scintillation material
- Cable malfunction

5.2 No Source Plateau

- Light leak, slide not sealed properly against true base
- Noisy PMT
- Noisy HV switch
- Poor PMT-to-scintillation light pipe interface

5.3 Excessive Background Count

- Light leak
- PMT malfunction
- Cable malfunction
- Noisy HV switch
- Instrument contaminated

**Model 43-10-6 Alpha-Beta Sample Counter
December 2015**

6. PARTS LIST

Ref. No.	Description	Part No.	Ref. No.	Description	Part No.
Model 43-10-6 Alpha Sample Counter			Switch Filter Board, Drawing 142 x 58		
UNIT	Assembled Model 43-10-6 Alpha/Beta Sample Counter	47-2794	BOARD	Assembled Switch Filter	5142-103
Assembly View, Drawing 142 x 113			CAPACITORS		
*	2 inch PM Tube Assembly	01-5919	C1-C2	0.0056μF, 3 KV	04-5522
*	Plastic Scintillator Material (EJ444L-2.20 x .010 ZnS)	01-5698	C3	0.0015μF, 3 KV	04-5518
*	Tube and Base Assembly	2142-128	RESISTORS		
*	M43-10 Connector Cap	7142-014	R1-R2	1 M, 1/4W, 5%	10-7028
*	M43-10-6 Sample Tray	7142-126	MISCELLANEOUS		
*	M43-10 Acrylic Disc	7142-002-01	*	CLVRLF-011-6809-000-599	18-8771
*	M43-10 Spacer Strip	7142-002-03	*	CONTACT-#1434	18-9124
*	M43-10 Adapter Plate	7142-003-01	Voltage Divider Board, Drawing 435 x 964		
*	M43-10 Case Top	7142-004-03	BOARD	Assembled Voltage Divider	5435-401
*	M43-10 Case Bottom	7142-004-04	CAPACITORS		
*	M43-10 Cap Gasket	7142-017	C1-C3	0.01μF, 2kV, C	04-5722
*	M43-10/43-17 Base Plate	7142-018	RESISTORS		
*	M43-10/43-17 Shaft	7142-019	R1-R12	4.75 MEG, 1/8W, 1%	10-7995
*	M43-10/43-17 Lifter	7142-020			
*	M43-10/43-17 Pin	7142-021			
*	2 inch X-TAL Foam Pad	7260-001-05			
*	Metallized polyester	01-5143			
*	M43-10-6 Filler Plate	7142-129			
*	M43-10 Bracket	7142-004-01			
*	M43-10 Cap	7142-004-02			
*	SWITCH-BZ-2RD-A2				
*	Micro Switch	08-6538			
*	KNOB-90 4 2G Pointer	08-6608			
*	RECEPT-UG568/U, C	4478-011			
*	BUMPER-Pad	21-9376			
*	Silicone Cord	22-9863			

**Model 43-10-6 Alpha-Beta Sample Counter
December 2015**

DRAWINGS AND DIAGRAMS

Assembly View, Drawing 142 x 113


Switch Filter Board, Drawing 142 x 58

Switch Filter Board Component Layout, Drawing 142 x 59

Voltage Divider Board, Drawing 435 x 964


Voltage Divider Board Component Layout, Drawing 435 x 965

REV #	ALTERATIONS	DATE	BY
2	CHANGED TUBE/SKT ASSY ECF-802	6-29-98	DDW
3	ECF # 1072	10/31/01	DSW
4	ECF # 2165	3/6/07	CMC
5	ECF # 2165	4/6/07	CMC




DWN	DATE	CHK	DATE	APP	DATE
CMC	4-6-07			<i>ROS GARD 7</i>	
PART NUM: 4142-181				SCALE: FULL <input type="checkbox"/> OTHER <input type="checkbox"/>	
TITLE M 43-10-6, M 43-10-6-1					
LUDLUM MEASUREMENTS, INC. 501 DAK STREET SWEETWATER, TEXAS 75556			SERIES	SHEET	
			142	113	


REVISIONS						
EFF	AUTHORITY	ZONE	LTR	DESCRIPTION	DATE	APPROVED




UPDATED	-	LUDLUM MEASUREMENTS INC.			
DR PH	10/20/92	TITLE: SWITCH			
CHK	CKB 07-JAN-91	FILTER BOARD			
DSCN PH	10/20/92	BOARD# 5142-103			
APPD	CS 11-6-01	SIZE	MODEL	SERIES	SHEET
NEXT HIGHER ASSY.	-	C	43-10	142	59
08:36:53	27-Jan-93	99142103	SHEET 1 OF 1		



<input checked="" type="checkbox"/> LUDLUM MEASUREMENTS INC. SWEETWATER, TX.			
DR	PW	10/20/92	TITLE: SWITCH FILTER BOARD
CHK	CKB	27-28-99	BOARD# 5142-103 BS142103
DSCN	PW	10/20/92	MODEL SERIES SHEET
APP	CKC	1-27-99	43-10 142 59
07:28:59	27-Jan-99	COMP SIDE <input type="checkbox"/>	SLDR SIDE <input type="checkbox"/> OUTLINE <input checked="" type="checkbox"/>
COMP PASTE <input type="checkbox"/>		COMP MASK <input type="checkbox"/> SLDR PASTE <input type="checkbox"/> SLDR MASK <input type="checkbox"/>	



		PO Box 810 501 Oak Street Sweetwater, Texas 79556 U.S.A. 1-800-622-0828	
Drawn: AC	05/07/2012	Title: VOLTAGE DIVIDER	
Design: RSS	05/07/2012	Model: VARIOUS	
		Board#: 5435-401	
Approve: <i>J.W.</i>	<i>10-22-12</i>	Sheet: 1 of 1	Series
Print Date: 10/10/2012 9:37:50 AM		Rev: 2	Sheet
W:\Projects\LMI\VoltageDividers\5435-401\Rev2\435401R2P1.SchDoc		435	964



LUDLUM
MEASUREMENTS, INC.

PO Box 810
501 Oak Street
Sweetwater, TX 79556
U.S.A. 1-800-622-0828

Title: VOLTAGE DIVIDER				
Drawn: AC	05/07/2012	Model: VARIOUS		
Design: RSS	05/07/2012	Board#: 5435-401		
Approve: <i>AW</i>	<i>10-22-12</i>	Rev: 2		
Print Date:		SCALE: 1.00	Series	Sheet
10/10/2012	9:37:52 AM	Top Overlay	435	965
W:\Projects\LMIVoltageDividers\5435-401\Rev2\435401R2_Manual.PcbDoc				