

# Model L-745

## X-ray Generator Output Detector



Ludlum Measurements, Inc.



### Introduction

The Ludlum Model L-745 X-ray Generator Output Detector provides a simple way to confirm X-ray generator performance. By connecting the detector directly to a digital scope, a storage scope, or a standard oscilloscope camera, the device will display intensity and time relationships of the X-ray beam to illustrate real-time operation. The detector supplies a minimum 200 mV signal at typical diagnostic settings (80 kVp /100 mA), which together with the better than 1 microsecond rise time, provides a true X-ray output pulse.

Using the device is straightforward. The detector, while connected to a standard oscilloscope, is simply placed in the X-ray beam. Examining the resulting waveform will enable a variety of potential X-ray issues to be identified. The general applications include measurement of exposure time, rectifier problems, cable or connector arcing, and general verification of X-ray operation. Model L-745 applications include documentation, verification, and troubleshooting of X-ray system operation. The detector is weighted for easy placement and uses a BNC-type connector.

### Specifications

Part Number 99-9214

**POWER SOURCE:** none required

**RISE TIME:** less than 1 microsecond

200 mV peak minimum

X-ray at 80 kVp and 100 mA with head 15 cm (6 in.) above the detector

**MINIMUM OSCILLOSCOPE INPUT IMPEDANCE:** 100,000 Ohms

**OSCILLOSCOPE IMPEDANCE:** typically 1 or 10 mOhm

**CONNECTOR:** BNC

**SIZE:** 51 x 32 x 20 mm (2.0 x 1.25 x 0.8 in.)

**WEIGHT:** 85 g (3 oz.)

**Option:** Cable Model L-40-1008-20

Ludlum Measurements, Inc. P.O. Box 810, Sweetwater, Texas 79556

**Web:** <http://www.ludlums.com> **Tel:** 800-622-0828 / 325-235-5494 / **Fax:** 325-235-4672 / **Email:** [sales@ludlums.com](mailto:sales@ludlums.com)

Note: specifications subject to change without notification. We are not responsible for errors or omissions.

March 2019