

# Model 43-10-1

## Alpha Sample Counter Head



Ludlum Measurements, Inc.



### Specifications

Part Number: 47-1305

**INDICATED USE:** alpha/beta sample counting

**SUGGESTED INSTRUMENTS:** Model 2929, Model 2223, Model 2224, Model 3030E

**DETECTOR TYPE:** scintillator, ZnS (Ag)

**DETECTOR OPERATING VOLTAGE:** 500–1200 volts

**CONNECTOR:** series "C" (others available)

**SAMPLE SIZE (maximum):** 5.1 x 0.9 cm (2 x 0.4 in.) (Dia x L)

**SAMPLE HOLDER:** anodized aluminum tray with sample ring to allow for 2.5 cm (1 in.) or 5.1 cm (2 in.) diameter x 0.32 cm (0.13 in.) deep samples

**PHOTOMULTIPLIER TUBE:** 5.1 cm (2 in.)

**WINDOW:** 0.4 mg/cm<sup>2</sup>

**EFFICIENCY (4π):** 37% for <sup>239</sup>Pu, 5% for <sup>14</sup>C, 27% for <sup>99</sup>Tc, 32% for <sup>230</sup>Th, 39% for <sup>238</sup>U, 29% for <sup>137</sup>Cs, 26% for <sup>99</sup>Sr/<sup>90</sup>Y

**BACKGROUND:** 3 cpm or less alpha, 80 cpm or less for beta-gamma

**CONSTRUCTION:** aluminum housing with beige powder coating

**TEMPERATURE RANGE:** -20 to 50 °C (-4 to 122 °F), may be certified for operation from -40 to 65 °C (-40 to 150 °F)

**SIZE:** 23.6 x 11.4 x 23.6 cm (9.3 x 4.5 x 9.3 in.) (H x W x L)

**WEIGHT:** 1.9 kg (4.1 lb)

### Optional Accessories

Sample Planchets

SIZE: 5.1 x 0.32 cm (2 x 0.125 in.) (Dia x Thickness)

Aluminum: Part Number 7525-371

Stainless Steel: Part Number 7525-371-01

(minimum order quantity 500)



Ludlum Measurements, Inc. P.O. Box 810, Sweetwater, Texas 79556

Web: <http://www.ludlums.com> Tel: 800-622-0828 / 325-235-5494 / Fax: 325-235-4672 / Email: [sales@ludlums.com](mailto:sales@ludlums.com)

Note: specifications subject to change without notification. We are not responsible for errors or omissions.

July 2014